



VALIDATE® Calibration Verification Test Set

Safety Data Sheet

according to Regulation (EC) No. 453/2010 and 29 CFR § 1910.1200

Revision date: Apr 19, 2024

Version: 9.0

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Product form : Mixture
Name : VALIDATE® Calibration Verification Test Set
Product code : 705

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Use of the substance/preparation : Intended for *in vitro* diagnostic use in the quantitative determination of linearity, calibration verification and verification of reportable range in automated, semi-automated and manual instrument systems.

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

LGC Clinical Diagnostics, Inc.
221 US Route 1
Cumberland Foreside, ME 04110 - USA
800-377-9684; 207-892-1300 Hours: 8:30 AM to 5:00 PM, Eastern Time, Monday through Thursday. 8:30 AM to 4:00 PM on Friday.

1.4. Emergency telephone number

Emergency number : Contact the CHEMTREC® Emergency Call Center for assistance with transportation or hazardous materials emergencies (24 hours/day, 7 days/week). Refer to LGC Clinical Diagnostics, Inc. customer number 798037.
Within USA and Canada: 1-800-424-9300 or +1 703-527-3887 (collect calls accepted)

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Not classified

Full text of H-phrases: see section 16.

Classification according to 2012 29 CFR § 1910.1200 [OSHA GHS]

Not classified

Classification according to WHMIS (Controlled Products Regulations)

Not classified

Adverse physicochemical, human health and environmental effects

None under normal conditions. However, product does contain material of human origin and should be considered potentially infectious.

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP) : No symbol required
Signal word (CLP) : Not Applicable
Hazard statements (CLP) : Not Applicable
Precautionary statements (CLP) : Not Applicable
EUH phrases : Not Applicable

Labelling according to 2012 29 CFR § 1910.1200 [OSHA]

Hazard pictograms (OSHA) : No symbol required
Signal word (OSHA) : Not Applicable
Hazard statements (OSHA) : Not Applicable
Precautionary statements (OSHA) : Not Applicable

2.3. Other hazards

None

PBT: not relevant – no registration required

vPvB: not relevant – no registration required

SECTION 3: Composition/information on ingredients

3.1. Substances

Not applicable

3.2. Mixtures

Name of Component	Product identifier	% in Mixture (see 2.1 for hazards for mixture)	Classification of Component according to Regulation (EC) No. 1272/2008 [CLP]
Human Serum Albumin	(CAS No.) 70024-90-7 EINECS:274-272-6	0.2-15	

Additional information: Contains human sourced and/or potentially infectious material.

VALIDATE® Calibration Verification Test Set

Safety Data Sheet

according to Regulation (EC) No. 453/2010 and 29 CFR § 1910.1200

SECTION 4: First aid measures

4.1. Description of first aid measures

- | | |
|---------------------------------------|---|
| First-aid measures after inhalation | : If symptoms are experienced remove source of contamination or move victim to fresh air. Obtain medical advice. |
| First-aid measures after skin contact | : No health effects expected. If irritation does occur, flush with lukewarm, gently flowing water for 5 minutes or until chemical is removed. |
| First-aid measures after eye contact | : No effects expected. If irritation does occur, flush contaminated eye(s) with lukewarm, gently flowing water for 5 minutes or until the chemical is removed. If irritation persists, obtain medical advice. |
| First-aid measures after ingestion | : If irritation or discomfort occurs, obtain medical advice immediately. |

4.2. Most important symptoms and effects, both acute and delayed

- | | |
|--------------------------------------|---|
| Symptoms/injuries after inhalation | : Overexposure to mists or spray may produce slight irritation of the respiratory system. |
| Symptoms/injuries after skin contact | : May cause slight irritation to the skin. |
| Symptoms/injuries after eye contact | : May cause slight irritation to eyes. |
| Symptoms/injuries after ingestion | : May cause gastric irritation. |

4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

SECTION 5: Firefighting measures

5.1. Extinguishing media

- | | |
|--------------------------------|--|
| Suitable extinguishing media: | : Alcohol-resistant foam and carbon dioxide. Extinguishing materials should be selected according to the surrounding area. |
| Unsuitable extinguishing media | : None known. |

5.2. Special hazards arising from the substance or mixture

- | | |
|-------------|--|
| Fire hazard | : No particular fire or explosion hazard. |
| Reactivity | : No dangerous reactions known under normal conditions of use. |

5.3. Advice for firefighters

- | | |
|---------------------------------------|---|
| Protective equipment for firefighters | : In case of fire: Wear self-contained breathing apparatus. Wear recommended personal protective equipment. |
|---------------------------------------|---|

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

6.1.1. For non-emergency personnel

- | | |
|----------------------|--|
| Protective equipment | : Wear suitable protective clothing, gloves and eye/face protection. |
|----------------------|--|

6.1.2. For emergency responders

- | | |
|----------------------|--|
| Protective equipment | : Wear suitable protective clothing, gloves and eye/face protection. |
| Emergency procedures | : Avoid contact with eyes and skin. Prevent product from entering drains. Prevent entry to sewers and public waters. |

6.2. Environmental precautions

Prevent product from entering drains.

6.3. Methods and material for containment and cleaning up

- | | |
|-------------------------|--|
| For containment | : Soak up spills with paper and absorbents. Wash with plenty of water and detergent. |
| Methods for cleaning up | : Soak up spills with paper and absorbents. Wash with plenty of water and detergent. Collect all waste in suitable and labelled containers and dispose according to local legislation. |

6.4. Reference to other sections

Refer to sections 8 and 13.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

- | | |
|-------------------------------|--|
| Precautions for safe handling | : Avoid all unnecessary exposure. Handle as potentially infectious material. |
| Hygiene measures | : Handle in accordance with good industrial hygiene and safety practice. Wash hands thoroughly after handling. Do not eat, drink or smoke when using this product. |

7.2. Conditions for safe storage, including any incompatibilities

- | | |
|------------------------|--|
| Technical measures: | : Provide local exhaust or general room ventilation. |
| Storage condition(s) | : Keep container tightly closed and dry. Store in well ventilated area. Keep cool. |
| Incompatible materials | : None known |

7.3. Specific end use(s)

Intended for in vitro diagnostic use in the quantitative determination of linearity, calibration verification and verification of reportable range in automated, semi-automated and manual instrument systems.

VALIDATE® Calibration Verification Test Set

Safety Data Sheet

according to Regulation (EC) No. 453/2010 and 29 CFR § 1910.1200

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

No additional information available

8.2. Exposure controls

- | | |
|----------------------------------|---|
| Appropriate engineering controls | : Either local exhaust or general room ventilation is usually required. |
| Personal protective equipment | : Gloves. Safety glasses. |



- | | |
|---------------------------------|--|
| Hand protection | : Wear suitable gloves resistant to chemical penetration. Consult your PPE provider for breakthrough times and the specific glove that is best for you based on your use patterns. |
| Eye protection | : As a minimum, safety glasses with side shields should be worn when handling this material. |
| Skin and Body protection | : If this material may come into contact with the body during handling and use, we recommend wearing appropriate protective clothing to prevent contact with the skin. |
| Respiratory protection | : With correct and proper use, and under normal conditions, breathing protection is not required. |
| Environmental exposure controls | : None known. |

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

- | | |
|--|---------------------------------|
| Physical state | : Liquid |
| Colour | : Colourless to slightly yellow |
| odour | : No data available |
| Odour threshold | : No data available |
| pH | : No data available |
| Melting point | : No data available |
| Solidification point | : No data available |
| Boiling point | : No data available |
| Flash point | : No data available |
| Relative evaporation rate (butylacetate=1) | : No data available |
| Flammability (solid, gas) | : Not applicable |
| Explosive limits | : No data available |
| Vapour pressure | : No data available |
| Relative vapour density at 20 °C | : No data available |
| Relative density | : No data available |
| Solubility | : Water: Soluble in water |
| Log Pow | : No data available |
| Log Kow | : No data available |
| Self ignition temperature | : No data available |
| Decomposition temperature | : No data available |
| Viscosity, kinematic | : No data available |
| Viscosity, dynamic | : No data available |
| Explosive properties | : Not explosive |
| Oxidising properties | : Not oxidizing |

9.2. Other information

- | | |
|----------------------------------|---------------------|
| Upper flammable Limit | : No data available |
| Lower Flammable Limit | : No data available |
| Sensitivity to mechanical impact | : No data available |
| Sensitivity to static discharge | : No data available |

VALIDATE® Calibration Verification Test Set

Safety Data Sheet

according to Regulation (EC) No. 453/2010 and 29 CFR § 1910.1200

SECTION 10: Stability and reactivity

10.1. Reactivity

No dangerous reactions known under normal conditions of use.

10.2. Chemical stability

Stable in use and storage conditions as recommended in Section 7.

10.3. Possibility of hazardous reactions

None under normal conditions.

10.4. Conditions to avoid

None under normal conditions.

10.5. Incompatible materials

None under normal conditions.

10.6. Hazardous decomposition products

No hazardous decomposition products known at room temperature.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Routes of Entry : Skin, Eyes, Inhalation and Ingestion
Acute toxicity : Not classified (Mixture not tested. Based on available data, the classification criteria are not met)

Human Serum Albumin CAS#70024-90-7

No data available

Skin corrosion/irritation : Not classified (Mixture not tested. Based on available data, the classification criteria are not met)
Serious eye damage/irritation : Not classified (Mixture not tested. Based on available data, the classification criteria are not met)
Respiratory or skin sensitisation : Not classified (Mixture not tested. Based on available data, the classification criteria are not met)
Germ cell mutagenicity : Not classified (Mixture not tested. Based on available data, the classification criteria are not met)
Carcinogenicity : Not classified (Mixture not tested. Based on available data, the classification criteria are not met).
No ingredient present at greater than 0.1% is listed in National Toxicology Program (NTP) Report on Carcinogens (12th Edition), the International Agency for Research on Cancer (IARC) Monographs (Volumes 1–105), or by OSHA.
Reproductive toxicity : Not classified (Mixture not tested. Based on available data, the classification criteria are not met)
Teratogenicity/Embryotoxicity : Not classified (Mixture not tested. Based on available data, the classification criteria are not met)
Specific target organ toxicity (single exposure) : Not classified (Mixture not tested. Based on available data, the classification criteria are not met)
Specific target organ toxicity (repeated exposure) : Not classified (Mixture not tested. Based on available data, the classification criteria are not met)
Aspiration hazard : Not classified (Not applicable)
Potential Adverse human health effects and symptoms : May cause slight temporary irritation.
Toxicologically Synergistic Materials : None known

SECTION 12: Ecological information

12.1. Toxicity

No additional information available

12.2. Persistence and degradability

No additional information available

12.3. Bioaccumulative potential

No additional information available

12.4. Mobility in soil

No additional information available

12.5. Results of PBT and vPvB assessment

VALIDATE® UC1 Calibration Verification Test Set

PBT: not relevant – no registration required

VALIDATE® Calibration Verification Test Set

Safety Data Sheet

according to Regulation (EC) No. 453/2010 and 29 CFR § 1910.1200

VALIDATE® UC1 Calibration Verification Test Set

vPvB: not relevant – no registration required

12.6. Other adverse effects

No additional information available

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste disposal recommendations : Dispose in a safe manner in accordance with local/national regulations. Dispose of at a licensed waste collection centre. Dispose of this material and its container at a hazardous or special waste collection point.

SECTION 14: Transport information

Not regulated for transport.

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1. EU-Regulations

Contains no REACH candidate substance

15.1.2. National regulations

Regional legislation : Canada: This product has been classified in accordance with the hazard criteria of the *Controlled Products Regulations* and the MSDS contains all the information required by the *Controlled Products Regulations*.
CAS#70024-90-7 is not listed on Canada's DSL List or Ingredient Disclosure List
USA: This product has been classified in accordance with the 2012 hazard criteria of the *Occupational Safety and Health Administration's (OSHA) Hazard Communication Standard (HCS)* and the SDS contains all the information required by the 29 CFR § 1910.1200.

15.2. Chemical safety assessment

CSA has not been established

SECTION 16: Other information

Sources of Key data : Data arise from reference works and literature.
Abbreviations and acronyms : CFR - Code of Federal Regulations
CLP - Classification, Labelling and Packaging
CSA – Chemical Safety Assessment
DPD – Dangerous Products Directive
EC - European Community
EC No. - European Commission number
EEC - European Economic Community
EN - European Norm
EU – European Union
GHS – Globally Harmonized System
h – hours
H-phrases – Hazard phrases
Log Pow – Log of the Octanol/Water Partition Coefficient
MSDS – Material Safety Data Sheet
No. – Number
OSHA - Occupational Safety and Health Administration
PBT - Persistent, Bioaccumulative and Toxic substance
PPE – Personal Protective Equipment
REACH - Registration, Evaluation, Authorisation and Restriction of Chemicals
R-phrases – Risk phrases
SDS – Safety Data Sheet
vPvB - Very Persistent and Very Bioaccumulative
WHMIS – Workplace Hazardous Material Information System

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product.